
KG-501

产品编号: D62000

CAS: 18228-17-6

分子式: C₁₇H₁₃ClNO₅P

纯度: 97%

InChi Key: InChIKey=RQAQWBFHPMSXKR-UHFFFAOYSA-N

Smiles: OP(O)(=O)Oc1cc2ccccc2cc1C(=O)Nc1ccc(Cl)cc1

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Epigenetic Reader Domain

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: KG-501 is a CREB inhibitor, with an IC₅₀ of 6.89 μM.