
YKL-06-061

产品编号: D62003

CAS: 2172617-15-9

分子式: C₃₀H₃₇N₇O₂

纯度: 98%

InChi Key: InChIKey=GVBAXIVNAHMIGH-UHFFFAOYSA-N

Smiles: COc1cc(ccc1Nc1ncc2CN(C(=O)N(C3CCC3)c2n1)c1c(C)cccc1C)N1CCN(C)CC1

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: SIK

溶解性: Soluble in DMSO

保存条件: 2-8°C

产品介绍: YKL-06-061 is a potent, selective, second-generation salt-inducible kinase (SIK) inhibitor with IC₅₀ values of 6.56 nM/1.77 nM/20.5 nM for SIK1/2/3, respectively