
(E/Z)-TG003

产品编号: D62006

CAS: 300801-52-9

分子式: C₁₃H₁₅NO₂S

纯度: 98%

InChi Key: InChIKey=BGVLELSCIHASRV-UHFFFAOYSA-N

Smiles: O=C(C=C1SC2=CC=C(OC)C=C2N1CC)C

外观: 固体粉末

作用通路: Apoptosis

作用靶点: BCL; CDK

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: (E/Z)-TG003 is a potent and ATP-competitive Cdc2-like kinase (Clk) inhibitor.