
BMS582949

产品编号: D62022

CAS: 623152-17-0

分子式: C₂₂H₂₆N₆O₂

纯度: 98%

InChi Key: InChIKey=GDTQLZHDDRBE-B-UHFFFAOYSA-N

Smiles: CCCNC(=O)c1cn2ncnc(Nc3cc(ccc3C)C(=O)NC3CC3)c2c1C

外观: 固体粉末

作用通路: MAPK

作用靶点: p38 MAPK

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: BMS-582949 is a potent and selective p38 mitogen-activated protein kinase (p38 MAPK) inhibitor with IC₅₀ of 13 nM, inhibiting both p38 kinase activity and activation of p38.