
Bay K 8644

产品编号: D62030

CAS: 71145-03-4

分子式: C₁₆H₁₅N₂O₄F₃

纯度: 98%

InChi Key: InChIKey=ZFLWDHHVRRZMEI-UHFFFAOYSA-N

Smiles: COC(=O)C1=C(C)NC(C)=C(C1C1=C(C=CC=C1)C(F)(F)F)[N+](=O)[O-]

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Calcium Channel

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Bay K 8644 (SQ 28,873) is a potent, selective activator of L-type Ca²⁺ channel with IC₅₀ of 17.3 nM.