

---

## KKL-35

产品编号: D62035

CAS: 865285-29-6

分子式: C<sub>15</sub>H<sub>9</sub>ClFN<sub>3</sub>O<sub>2</sub>

纯度: 98%

InChi Key: InChIKey=ZIICPNCCHIUJSX-UHFFFAOYSA-N

Smiles: Fc1ccc(cc1)-c1nnc(NC(=O)c2ccc(Cl)cc2)o1

外观: 固体粉末

作用通路: Microbiology/Virology

作用靶点: Antibacterial

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: KKL-35 is a trans-translation tagging reaction inhibitor with an IC<sub>50</sub> of 0.9 μM