
DLinDMA

产品编号: D62039

CAS: 871258-12-7

分子式: C₄₁H₇₇NO₂

纯度: 98%

InChi Key: InChiKey=NFQBIAXADRDUGK-KWXKLSQISA-N

Smiles: CCCCC/C=C\C/C=C\CCCCCCCCOCC(OCCCCCCCC/C=C\C/C=C\CCCC)CN(C)C

外观: 固体粉末

作用靶点: Others

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: DLinDMA, a ionizable cationic lipid, is a key lipid component of stable nucleic acid lipid particles (SNALPs) as a benchmark. DLinDMA is used for siRNA delivery