
SF2523

产品编号: D62062

CAS: 1174428-47-7

分子式: C₁₉H₁₇NO₅S

纯度: 98%

InChi Key: InChIKey=BYTKNUOMWLJVNQ-UHFFFAOYSA-N

Smiles: O=c1cc(N2CCOCC2)oc2c1scc2c1ccc2OCCOc2c1

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: DNA-PK; Epigenetic Reader Domain; mTOR; PI3K

溶解性: Soluble in DMSO、Ethanol、H₂O

保存条件: -20°C

产品介绍: SF2523 is a highly selective and potent inhibitor.