
S 38093 HCl

产品编号: D62069

CAS: 1222097-72-4

分子式: C₁₇H₂₅ClN₂O₂

纯度: 98%

InChi Key: InChIKey=AFSVOZDCVFWFG-UHFFFAOYSA-N

Smiles: Cl.NC(=O)c1ccc(cc1)OCCCN2CC3CCCC3C2

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Histamine Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: S 38093 is a novel brain-penetrant antagonist/inverse agonist of H3 receptors