
BQR-695

产品编号: D62093

CAS: 1513879-21-4

分子式: C₁₉H₂₀N₄O₃

纯度: 98%

InChi Key: InChIKey=LYPECULYCGFOIDA-UHFFFAOYSA-N

Smiles: CNC(=O)CNc1cnc2ccc(cc2n1)c1cc(c(cc1)OC)OC

外观: 固体粉末

作用通路: Microbiology/Virology

作用靶点: Parasite; PI3K; PI4K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: BQR-695 is a PI4KIIIβ inhibitor with IC₅₀s of 80 and 3.5 nM for human PI4KIIIβ and Plasmodium variant of PI4KIIIβ, respectively.