
TP0463518

产品编号: D62094

CAS: 1558021-37-6

分子式: C₂₀H₁₈ClN₃O₆

纯度: 98%

InChi Key: InChIKey=HMMHKGLPKAQOOH-UHFFFAOYSA-N

Smiles: OC(=O)CNC(=O)C1=C(O)CCN(Cc2ccc(Oc3ccc(Cl)cc3)nc2)C1=O

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: HIF/HIF Prolyl-Hydroxylase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: TP-0463518 is a highly potent HIF prolyl hydroxylase (PHD) inhibitor (IC₅₀s: 13 nM and 18 nM for human and rat PHD2, respectively).