
P110δ-IN-1

产品编号: D62096

CAS: 1595129-71-7

分子式: C₃₁H₄₀N₈O₃S

纯度: 98%

InChi Key: InChIKey=PGTXVGREXBMCCY-UHFFFAOYSA-N

Smiles: O=S(C1=NC2=CC=CC=C2N1C3=NC(N4CCOCC4)=NC(NC(C)(C)CC5=CC=CC=C5C6CCN(C)CC6)=N3)(C)=O

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: P110δ-IN-1 (PWT-143) is an effective and selective inhibitor of P110δ (IC₅₀: 8.4 nM).