
ARN-3236

产品编号: D62101

CAS: 1613710-01-2

分子式: C₁₉H₁₆N₂O₂S

纯度: 98%

InChi Key: InChIKey=WEHOIIGXTMKVRG-UHFFFAOYSA-N

Smiles: COc1ccc(-c2c[nH]c3nccc(-c4ccsc4)c23)c(OC)c1

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: SIK

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ARN-3236 is an oral active and selective inhibitor of salt-inducible kinase 2 (SIK2), with IC₅₀s of <1 nM, 21.63 nM and 6.63 nM for SIK2, SIK1 and SIK3, respectively. Has anti-cancer activity