
GLPG1837

产品编号: D62110

CAS: 1654725-02-6

分子式: C₁₆H₂₀N₄O₃S

纯度: 98%

InChi Key: InChIKey=GHTGYZMBQPXTCQ-UHFFFAOYSA-N

Smiles: CC1(C)Cc2c(sc(NC(=O)c3cc[nH]n3)c2C(N)=O)C(C)(C)O1

外观: 固体粉末

作用通路: Autophagy

作用靶点: Autophagy; CFTR

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: GLPG1837 is an effective CFTR potentiator, with EC₅₀s of 3 nM and 339 nM for F508del and G551D CFTR, respectively.