
TAS6417

产品编号: D62111

CAS: 1661854-97-2

分子式: C₂₃H₂₀N₆O

纯度: 98%

InChi Key: InChIKey=MKCYPWYURWOKST-INIZCTEOSA-N

Smiles: CC1=C[C@@H](CN2C1=C(C3=C(N=CN=C32)N)C4=CC5=CC=CC=C5N=C4)NC(=O)C=C

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: Apoptosis; EGFR

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: TAS6417 is an EGFR inhibitor and is an efficacious drug candidate for patients with NSCLC (IC₅₀: ranging from 1.1-8.0 nM).