
THZ531

产品编号: D62116

CAS: 1702809-17-3

分子式: C₃₀H₃₂ClN₇O₂

纯度: 98%

InChi Key: InChIKey=RUBYHLPRZRM TJ O-MOVYNIQHSA-N

Smiles: O=C(NC1=CC=C(C(N2C[C@H](NC3=NC=C(Cl)C(C4=CNC5=C4C=CC=C5)=N3)CCC2)=O)C=C1)/C=C/CN(C)C

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: CDK

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: THZ531 is a selective and covalent inhibitor of both CDK12 and CDK13 with IC₅₀s of 158 nM and 69 nM, respectively