
SR-4370

产品编号: D62124

CAS: 1816294-67-3

分子式: C₁₇H₁₈F₂N₂O

纯度: 98%

InChi Key: InChIKey=OSQKWTZHYXRTBG-UHFFFAOYSA-N

Smiles: Fc1c(cccc1F)c2ccc(cc2)C(=O)NNCCCC

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: HDAC

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: SR-4370 is an inhibitor of HDAC, with IC₅₀s of 0.13 μM, 0.58 μM, 0.006 μM, 2.3 μM, and 3.4 μM for HDAC1, HDAC2, HDAC3, HDAC8, and HDAC6, respectively.