
NCT-503

产品编号: D62131

CAS: 1916571-90-8

分子式: C₂₀H₂₃F₃N₄S

纯度: 98%

InChi Key: InChIKey=PJNSZIQUFLWRLH-UHFFFAOYSA-N

Smiles: Cc1cc(nc(c1)NC(=S)N1CCN(CC1)Cc1ccc(cc1)C(F)(F)F)C

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: Dehydrogenase; FLT; TAM Receptor

溶解性: Soluble in DMSO、Ethanol

保存条件: 2-8°C

产品介绍: NCT-503 is a phosphoglycerate dehydrogenase (PHGDH) inhibitor with an IC₅₀ of 2.5 μM.