
WM-8014

产品编号: D62137

CAS: 2055397-18-5

分子式: C₂₀H₁₇FN₂O₃S

纯度: 98%

InChi Key: InChIKey=PHHZKBMVRPULAW-UHFFFAOYSA-N

Smiles: O=S(=O)(NNC(=O)c1cc(cc(C)c1F)c2ccccc2)c3ccccc3

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Acetyltransferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: MOZ-IN-3 is an inhibitor of MOZ(IC₅₀=55 nM), a member of histone acetyltransferases.