
WM-1119

产品编号: D62138

CAS: 2055397-28-7

分子式: C₁₈H₁₃F₂N₃O₃S

纯度: 98%

InChi Key: InChiKey=QLXULUNLCRKWRD-UHFFFAOYSA-N

Smiles: Fc1cccc1S(=O)(=O)NNC(=O)c2cc(cc(F)c2)c3ccccc3

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Acetyltransferase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: WM-1119 is a highly potent and selective KAT6A inhibitor, with an IC₅₀ of 0.25 μM for KAT6A in lymphoma cells, the binding K_D values of WM-1119 with KAT6A, KAT5 and KAT7 are 2 nM, 2.2 μM, 0.5 μM , respectively