
WNK-IN-11

产品编号: D62144

CAS: 2123489-30-3

分子式: C₂₁H₂₁Cl₂N₅O₅

纯度: 98%

InChi Key: InChIKey=MVXAYIXYYOVALX-UHFFFAOYSA-N

Smiles: CNc1nc(cs1)-c1cc(C(=O)N2CCN(Cc3ccc(Cl)cc3)CC2)c(Cl)cn1

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: Serine/threonin kinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: WNK-IN-11 is an allosteric With-No-Lysine (WNK) kinase inhibitor, with an IC₅₀ of 4 nM for WNK1.