
SGC-GAK-1

产品编号: D62154

CAS: 2226517-76-4

分子式: C₁₈H₁₇BrN₂O₃

纯度: 98%

InChi Key: InChIKey=AUOSKLDNVNGKRR-UHFFFAOYSA-N

Smiles: COc1cc(Nc2ccnc3ccc(Br)cc23)cc(OC)c1OC

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: Serine/threonin kinase

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: SGC-GAK-1 is a potent, selective cyclin G-associated kinase (GAK) inhibitor with a K_i of 3.1 nM. SGC-GAK-1 is a chemical probe for GAK