
Tetracycline

产品编号: D21042

CAS: 60-54-8

保存条件: Store at -20°C for 3 year(Powder); In DMSO or others solvent store at 2-4°C for two weeks, at -20°C for six months.

产品介绍: CAS: 60-54-8

分子式: C₂₂H₂₄N₂O₈

分子量: 444.43

MDL: MFCD00151232

纯度: ≥90% (HPLC)

作用靶点: Bacterial; Antibiotic;

作用通路: Anti-infection;

产品描述: 四环素 (Tetracycline) 是一种广谱抗生素, 对多种革兰氏阳性和革兰氏阴性细菌有抑制活性。

PRODUCT SPECIFIC PUBLICATIONS

[IF=8.5] Ying Zhou. et al. Efficient fluorescence and electrochemiluminescence dual-signal au nanoclusters-based portable antibiotic testing platform with super-wide detection range. FOOD CHEM. 2025 Sep;485:144339 ; . 40306061