
AZD4635

产品编号: D62239

CAS: 1321514-06-0

分子式: C₁₅H₁₁ClFN₅

纯度: 98%

InChi Key: InChIKey=NCWQLHHDGDXIJN-UHFFFAOYSA-N

Smiles: CC1=NC(=CC(=C1)C2=C(N=C(N=N2)N)C3=CC=C(C=C3)F)Cl

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Adenosine Receptor

溶解性: Soluble in DMSO、Ethanol

保存条件: -20°C

产品介绍: AZD4635 (HTL1071) is a new-type adenosine 2A receptor (A2AR) inhibitor(Ki=1.7 nM).