
HM03

产品编号: D62244

CAS: 500565-15-1

分子式: $C_{26}H_{27}ClN_4O_2$

纯度: 98%

InChi Key: InChIKey=SUSDGTMJKOGWSZ-UHFFFAOYSA-N

Smiles: COc1ccc2nc3cc(Cl)ccc3c(Nc3ccc(O)c(CN4CCN(C)CC4)c3)c2c1

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: HSP

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: HM03 is a potent and selective HSPA5 (Heat shock 70kDa protein 5, also known as Bip, Grp78) inhibitor. HM03 has anticancer activity