
ETC-206

产品编号: D62245

CAS: 1464151-33-4

分子式: $C_{25}H_{20}N_4O_2$

纯度: 98%

InChi Key: InChIKey=FWRFPHJSGLYXTD-UHFFFAOYSA-N

Smiles: O=C(N1CCOCC1)c1ccc(cc1)-c1ccc2ncc(-c3ccc(cc3)C#N)n2c1

外观: 固体粉末

作用通路: MAPK

作用靶点: MNK

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Tinodasertib (ETC-206) is a selective MNK1 and MNK2 inhibitor with IC50s of 64 nM and 86 nM, respectively.