

---

## LY255283

产品编号: D62253

CAS: 117690-79-6

分子式: C<sub>19</sub>H<sub>28</sub>N<sub>4</sub>O<sub>3</sub>

纯度: 98%

InChi Key: InChIKey=WCGXJPFHTHQNJL-UHFFFAOYSA-N

Smiles: CCc1cc(C(C)=O)c(O)cc1OCCCCC(C)(C)c1nn[nH]n1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Leukotriene Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: LY255283 is a LTB<sub>4</sub> receptor (BLT<sub>2</sub>) antagonist, with an IC<sub>50</sub> of ~100 nM for [3H]LTB<sub>4</sub> binding to guinea pig lung membranes