
ML-180

产品编号: D62266

CAS: 863588-32-3

分子式: C₂₀H₂₅ClN₄O₂

纯度: 98%

InChi Key: InChIKey=CHMQQIFPOIZPIJ-UHFFFAOYSA-N

Smiles: Clc1cccc(c1)N1CCN(CC1)c1cc(=O)n(C2CCCCC2)c(=O)[nH]1

外观: 固体粉末

作用靶点: Others

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ML-180 (SR1848) is a potent orphan nuclear receptor liver receptor homolog 1 (LRH-1; NR5A2) inverse agonist with an IC₅₀ of 3.7 μM. ML-180 is inactive for steroidogenic factor-1 (SF-1; NR5A1; IC₅₀>10 μM). ML-180 has the potential for LRH-1-dependent cancers