
TTK21

产品编号: D62272

CAS: 709676-56-2

分子式: C₁₇H₁₅ClF₃NO₂

纯度: 98%

InChi Key: InChIKey=AQJBXYBDNZHZRE-UHFFFAOYSA-N

Smiles: O=C(NC1=CC=C(Cl)C(=C1)C(F)(F)F)C=2C=CC=CC2OCCC

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Epigenetic Reader Domain

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: TTK21 is an activator of CBP/p300 histone acetyltransferase activity.