
APS6-45

产品编号: D62293

CAS: 2188236-41-9

分子式: C₂₃H₁₆F₈N₄O₃

纯度: 98%

InChi Key: InChIKey=MVUPJRHPAMCBQJ-UHFFFAOYSA-N

Smiles: O=C(NC)C1=NC=CC(OC2=CC=C(NC(NC3=CC(C(C(F)(F)F)(F)C(F)(F)F)=CC=C3F)=O)C=C2)=C1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: MAPK; Ras

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: APS6-45 is an orally active tumor-calibrated inhibitor (TCI). APS6-45 inhibits RAS/MAPK signaling and exhibits antitumor activity