

---

## Eliapixant

产品编号: D62397

CAS: 1948229-21-7

分子式: C<sub>22</sub>H<sub>21</sub>F<sub>3</sub>N<sub>4</sub>O<sub>3</sub>S

纯度: 98%

InChi Key: InChIKey=GTQPEQGDLVSFJO-CXAGYDPISA-N

Smiles: C[C@H](c1cnc(C(F)(F)F)nc1)NC(c1cc(-c2ncc(C)s2)cc(O[C@H]2COCC2)c1)=O

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: P2X Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: Eliapixant (BAY 1817080) is a potent and selective antagonist of P2X<sub>3</sub> receptor, with an IC<sub>50</sub> of 8 nM. Eliapixant can be used for the research of refractory chronic cough