
5-HT3 antagonist-4i

产品编号: D62510

CAS: 930478-88-9

分子式: C₁₆H₁₂ClN₃O

纯度: 98%

InChi Key: InChIKey=SCIBFCFFUQXLCN-UHFFFAOYSA-N

Smiles: Cc1c(Cl)cccc1NC(=O)c1cnc2ccccc2n1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: 5-HT Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: 4i is a 5-HT3 receptor antagonist, It modulates the serotonergic system.