
PDGFR Tyrosine Kinase Inhibitor III

产品编号: D62532

CAS: 205254-94-0

分子式: C₂₇H₂₇N₅O₄

纯度: 98%

InChi Key: InChIKey=INTPTKHSGKBHHW-UHFFFAOYSA-N

Smiles: O=C(NC1=CC=C(OC=CC=CC2)C=C1)N3CCN(C=4N=CN=C5C=C(OC)C(OC)=CC54)CC3

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: PDGFR

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: PDGFR Tyrosine Kinase Inhibitor III (PDGF Receptor Tyrosine Kinase Inhibitor III), a multikinase inhibitor, inhibits PDGFR, EGFR, FGFR, PKA, and PKC, respectively. PDGFR Tyrosine Kinase Inhibitor III can be used for the research of amyotrophic lateral sclerosis