
TC-AQP1-1

产品编号: D62538

CAS: 37710-81-9

分子式: C₁₂H₁₀O₄

纯度: 98%

InChi Key: InChIKey=KRXUBZPHAPGHPE-UHFFFAOYSA-N

Smiles: O=C(O)C=CC1=CC=CC(C=CC(=O)O)=C1

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Aquaporin

保存条件: -20°C

产品介绍: TC AQP1 1 (m-Phenylenediacrylic acid) is a potent AQP1 inhibitor. TC AQP1 1 inhibits AQP1-mediated water flux in oocytes, with an IC₅₀ of 8 μM