
p38 α inhibitor 3

产品编号: D62544

CAS: 260428-69-1

分子式: C₁₉H₂₀FNO

纯度: 97%

InChi Key: InChiKey=VALGZDLHWXYOBZ-UHFFFAOYSA-N

Smiles: O=C(C1=CC=C(F)C=C1)N2CCC(CC=3C=CC=CC3)CC2

外观: 固体粉末

作用通路: MAPK

作用靶点: p38 MAPK

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: p38 α inhibitor 3 (Comp G7) is a p38 α inhibitor that blocks the effectiveness of myoblast differentiation