
SP-100030

产品编号: D62565

CAS: 154563-54-9

分子式: C₁₄H₅ClF₉N₃O

纯度: 98%

InChi Key: InChIKey=GHMKQBWHQMXSH-UHFFFAOYSA-N

Smiles: O=C(NC=1C=C(C=C(C1)C(F)(F)F)C(F)(F)F)C2=CN=C(Cl)N=C2C(F)(F)F

外观: 固体粉末

作用靶点: NF-κB

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: SP-100030 is a potent NF-κB and activator protein-1 (AP-1) double inhibitor (IC₅₀s=50 and 50 nM, respectively). SP-100030 inhibits IL-2, IL-8, and TNF-alpha production in Jurkat and other T cell lines. SP-100030 decreases murine collagen-induced arthritis (CIA).