
HIF1-IN-3

产品编号: D62589

CAS: 333314-79-7

分子式: C₂₆H₂₄N₂O₃

纯度: 98%

InChi Key: InChIKey=RFCCQJJHFAITGH-UHFFFAOYSA-N

Smiles: COc1ccc(C(c(ccc2cccnc22)c2O)NC(Cc2ccccc2)=O)cc1

外观: 固体粉末

作用通路: Angiogenesis

作用靶点: HIF

溶解性: Soluble in DMSO

保存条件: 2-8°C

产品介绍: HIF1-IN-3 (compound F4) is a potent HIF1 inhibitor with an EC₅₀ value of 0.9 μM. HIF1-IN-3 can be used for researching anticancer