
NF-κB activator 1

产品编号: D62611

CAS: 2387524-59-4

分子式: C₁₆H₁₁FN₂O₂S

纯度: 98%

InChi Key: InChIKey=SSUFECODNVILSX-UHFFFAOYSA-N

Smiles: O=C(O)C1=CC=CC(=C1)C=2N=C(SC2)NC=3C=CC=CC3F

外观: 固体粉末

作用靶点: NF-κB

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: NF-κB activator 1 can activate the nf-kappa B gene activator, EC₅₀ of 0.9 microns. NF-κB is widely used in eukaryotic cells as a gene regulating cell proliferation and cell survival. NF-κB activator 1 induces superoxide dismutase (SOD) 2.1 mRNA expression.