
ST-1006

产品编号: D62634

CAS: 1196994-11-2

分子式: C₁₆H₂₀Cl₂N₆

纯度: 98%

InChi Key: InChIKey=ICAQZJCLBRKWBG-UHFFFAOYSA-N

Smiles: ClC1=CC=CC(Cl)=C1CNC2=NC(=NC(=C2)N3CCN(C)CC3)N

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Histamine Receptor

保存条件: -20°C

产品介绍: ST-1006 is a potent histamine H₄ receptor agonist (pK_i: 7.94). ST-1006 has anti-inflammatory effects. ST-1006 is a potent inducer of basophil migration and induces basophil migration. ST-1006 is a potent inducer of basophil migration. ST-1006 has anti-inflammatory and anti-pruritic effects.