
BAY-4931

产品编号: D62636

CAS: 423150-91-8

分子式: C₂₂H₁₆ClN₃O₄

纯度: 98%

InChi Key: InChIKey=WFEBALYYDASKIV-UHFFFAOYSA-N

Smiles: O=C(NC1=CC=C2OC(=NC2=C1)C3=CC=C(C=C3)CC)C4=CC(=CC=C4Cl)N(=O)=O

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: PPAR

保存条件: -20°C

产品介绍: BAY-4931 is a powerful, covalent, and selective inverse-agonist of PPAR γ , exhibiting an IC₅₀ value of 0.17 nM.