
A2A receptor antagonist 1

产品编号: D62649

CAS: 443103-97-7

分子式: C₁₆H₁₂FN₅O

纯度: 98%

InChi Key: InChIKey=JEEJMSUHUZNTCD-UHFFFAOYSA-N

Smiles: NC1=NC(C2=CC=CO2)=C3C(N(CC4=CC=CC=C4F)N=C3)=N1

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Adenosine Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: A2A receptor antagonist 1 is an antagonist of both adenosine A2A receptor and A1 receptor with K_is of 4 and 264 nM, respectively.