
DPTN dihydrochloride

产品编号: D62655

CAS: 325767-87-1

分子式: C₂₂H₂₀Cl₂N₄O₂S

纯度: 98%

InChi Key: InChIKey=GTWKWBXTCBUKTR-UHFFFAOYSA-N

Smiles: O=C(C1=CC=CN=C1)NC2=NC(C3=CC(C)=CC(C)=C3)=C(C4=CC=NC=C4)S2.Cl.Cl

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: Adenosine Receptor

保存条件: -20°C

产品介绍: DPTN dihydrochloride is a potent and selective A₃AR antagonist with K_i values of 1.65, 9.61 and 8.53 nM in human, mouse and rat, respectively.