

**bsm-41832M**

**[ Primary Antibody ]**

## **MH (Coating) Mouse mAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG2b

**Applications:** ELISA

**Clonality:** Monoclonal

**CloneNo.:** 3C12

**Reactivity:** (predicted: Mycoplasma hominis)

**Target:** MH (Coating)

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.