

---

## MH (Coating) Mouse mAb

Catalog Number: bsm-41832M

Target Protein: MH (Coating)

Concentration: 1mg/ml

Form: Liquid

Host: Mouse

Clonality: Monoclonal

Clone No.: 3C12

Isotype: IgG2b

Applications: ELISA

Reactivity: (predicted:Mycoplasma hominis)

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.