

**bsm-41833M**

**[ Primary Antibody ]**

## **MH (Detection) Mouse mAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG2b

**Applications:** ELISA (1:5000-10000)

**Clonality:** Monoclonal

**CloneNo.:** 7H5

**Reactivity:** (predicted: Mycoplasma hominis)

**Target:** MH (Detection)

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles