

bsm-41833M

[Primary Antibody]

MH (Detection) Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG2b

Applications: ELISA

Clonality: Monoclonal

CloneNo.: 7H5

Reactivity: (predicted: Mycoplasma hominis)

Target: MH (Detection)

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles