

---

## Human CEA Mouse mAb, BF488 conjugated

Cat. No:	V1904-BF488
Name:	Human CEA Mouse mAb, BF488 conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C9
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted:&nbsp;Human)
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Antigen:	Recombinant Human CEACAM5 Protein