

---

## Human CEA Mouse mAb, BF488 conjugated

|                |   |
|----------------|---|
| Cat. No:       | V1904-BF488   |
| Name:          | Human CEA Mouse mAb, BF488 conjugated   |
| Storage:       | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:    | Monoclonal  |
| Clone No:      | 6C9   |
| Host/Source:   | Mouse   |
| Subtype:       | IgG1  |
| Reactivity:    | (predicted:Human)   |
| Concentration: | 1mg/ml  |
| Buffer:        | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.               |
| Antigen:       | Recombinant Human CEACAM5 Protein   |