
Human CEA Mouse mAb, AP conjugated

Cat. No:	V1904-AP
Name:	Human CEA Mouse mAb, AP conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6C9
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Buffer:	0.01M PBS (pH7.4).
Antigen:	Recombinant Human CEACAM5 Protein