

bsm-0365M-Cy5

[Secondary Antibodies]

human IgG-Fc, Cy5 conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Target: human IgG-Fc

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.

Applications: IF (1:100-1000)

Flow-Cyt (1:100-1000)

ICC/IF (1:100-1000)

Excitation Spectrum: 649nm

Emission spectrum: 670nm

Reactivity: Human

Subcellular Location: Secreted